

Chemical Tests for Basic Nutrients

General Instruction: Use *small* samples of food. A piece this big: **0** is enough. Grind the food well in the mortar to increase the surface area of the sample. If the food is a *liquid*, put $\frac{1}{2}$ **cm** (no more!) in a test tube. *Excessive amounts skew the results of the tests!*

Starch: Place a sample of food in a space on a *spot plate*. Place a drop of **iodine** on the food.

- * Treated spot turns black, blue, or purple
- + Color change after several minutes, or around edges only
- o No color change, or spot turns iodine-color (brown)

Fat: Rub or squash a small amount of the food onto a *paper towel*. Place the towel on the *hot plate*, set on low, and examine after fifteen minutes.

- * Large grease spot under the food
- + Small grease spot
- o No grease spot

Note: most foods will leave a water spot which initially may appear as a grease spot, but the water will evaporate and the grease will not.

Protein: Place a small amount of food in a test tube and add just enough **10% NaOH** (Sodium hydroxide) to cover it. Then add **0.5% CuSO₄** (copper sulfate) drop-by-drop and count the number of drops you add.

- * Liquid becomes purple or pink after 1 - 10 drops
- + Liquid becomes purple or pink after 10 - 20 drops
- o No color change

Maltose: (Sugars) Place a small amount of food in a test tube and add enough **Benedict's** solution to cover it. Place the tube in *boiling water bath* or over a *Bunsen burner*. Watch for color changes in the liquid (not in the food).

- ** Liquid turns brown or red
- * Liquid turns orange or yellow
- + Liquid turns green
- o No color change (liquid remains blue).

Glucose: Place ground food sample in a test tube. Add **1 cm of water** and shake the tube to dissolve or mix the food. Dip a **Clinistix** into the solution and match it to the label on the bottle.

- ** Matches dark purple square on the bottle
- * Matches medium lavender square on the bottle
- + Matches pale lavender square on the bottle
- o Matches dark pink square on the bottle

Name: _____

Water: Place a food sample in a *dry test tube*. Heat gently over a *Bunsen burner flame*. Look for drops of moisture forming on the inside of the test tube.

- * Moisture appears
- o No moisture

Vitamin C: This test can only be done with *liquid foods* or foods which can be completely dissolved. Place **1.0 ml of 0.1% Indophenol** in a test tube. Add your food solution drop-by-drop and count the number of drops you add. At first the Indophenol may turn pink, which does not count – keep adding until the Indophenol becomes colorless. If it is still pink or blue after 20 drops, stop.

- * Liquid turns colorless after 1 - 10 drops
- + Liquid turns colorless after 11 - 20 drops
- o No color change (liquid remains pink or blue).

Food Tests	Starch	Fat	Protein	Sugar	Glucose	Water	Vit C
1)							
2)							
3)							
4)							
5)							
Unknown #							

Questions

1. Did any of the foods tested contain all four basic nutrients? _____
2. If any, which _____
Did you test any foods which seemed to have no nutrients Value? _____
3. Which food seemed to be the least nutritious? _____
Is this food really “worthless”? _____ Explain: _____

If vitamin C was tested, what seems to be the effect of the following:
Passage of time (storage)? _____
Cooking? _____