

Btec 101 Batch Record Activity

The assignment for your group is to refer to the example batch record to create a batch record for the following formula to make LB media:

- Bacto-tryptone 10g
- Bacto-yeast extract 5g
- NaCl 10g
- QS to 1 liter

Each lab session, bring a printed draft copy of your group's LB batch record. Use it when making your media, signing off where appropriate. Those steps that were incomplete or need to be edited, should be corrected and printed out in a new draft batch record for the next lab period. Repeat until you have a final version to turn in.

Your group's Word processed batch record is due on March 1 (day class),
March 2 (night class).